

# TQC SHEEN SCRUB ABRASION AND WASHABILITY TESTER STANDARDS CONFIGURATION TABLE

x = mandatory  
o = optional

Art. No.	Description	DIN 53778	ASTM D2486	ASTM D4213-92	ASTM D4828	ISO 11998 / EN13300	ASTM D3450	ASTM D4213 (from 1996)	MEK test / DIN 55654 / ISO 105-X1254.2.1
<b>Scrub Test Basic Device</b>									
AB5000	Scrub, Abrasion and Washability Test 230 V	x	x	x	x	x	x	x	x
AB5005	Scrub, Abrasion and Washability Test 100 V	x	x	x	x	x	x	x	x
AB5210	Scrub, Abrasion and Washability Test Basic 230 V		x	x	x		x	x	x
AB5215	Scrub, Abrasion and Washability Test Basic 100 V		x	x	x		x	x	x
<b>Test beds</b>									
AB5196	Reservoir low, half	x	x	x	x	x	x	x	x
AB5176	Reservoir low, full	o	o	o	o	o	o	o	o
<b>Clamping frames</b>									
AB5045	Wide channel single frame for low reservoir	x	x			x			x
AB5197	Double narrow channel frame for low reservoir	o	o	x	x	o	x	x	o
<b>Tool adapter</b> (WARNING: The adapter must be compatible with the clamping frame)									
AB5041	Tool adapter "Double" for tools AB5012/5010/5011/5013/5020	x	x			x			x
AB5042	Tool adapter "Single" for tools AB5012/5010/5011/5013/5020	o	o			o			o
AB5043	Tool adapter "Single" for tool AB5050			x	x		x	x	
AB5044	Tool adapter "Double" for tool AB5055								
<b>Tools</b>									
AB5010	Wild Boar Brush (250g/8.82oz) DIN 53778	x							
AB5011	Nylon Brush (454g/16.01oz) ASTM D2486		x						
AB5012	Sponge (232g/8.12oz)								
AB5013	Abrasive pad (135g/4.76oz) ISO 11998					x			
AB5020	Universal holder								
AB5050	Sponge tool, large, conform ASTM D4828 / D3450 (500g)(excl weights)			x	x		x	x	
AB5055	Cloth holder tool for detergents test (norms)								
AB5060	MEK Test / Crock meter								x
<b>Accessories</b>									
AB5014	Weight for Washability test 60g								
AB5015	Weight for Washability test 100g								
AB5030	Weight for Washability test 1000g						o		
AB5040	Weight for Washability test 500gr						o		
AB5032	Tool stabilization kit for heavy weights	o	x	o	o	o	x	o	
AB5025	Metal shims for ASTM D2486, for reservoir half (set of 4)		x						
AB5027	Metal shims for ASTM D2486, for reservoir full (set of 4)		o						
AB5180	Glass base conform ASTM D2486 (set of 2)	x	x	x	x		x	x	

Art. No.	Description	DIN 53778	ASTM D2486	ASTM D4213-92	ASTM D4828	ISO 11998 / EN13300	ASTM D3450	ASTM D4213 (from 1996)	MEK test / DIN 55654 / ISO 105-X1254.2.1
----------	-------------	-----------	------------	---------------	------------	---------------------	------------	------------------------	--

**Consumables**

<b>AB5016</b>	Abrasive pads for tool AB5013, set 5 pcs (Weight each appr. 3 g.)					<b>x</b>			
<b>AB5016</b>	Abrasive pads for tool AB5050, set 5 pcs			<b>x</b>				<b>x</b>	
<b>AB5017</b>	Wild boar brush for tool AB5010	x							
<b>AB5018</b>	Nylon brush for tool AB5011		<b>x</b>						
<b>AB5018</b>	Sponge pads for tool AB5050, set 5 pcs			<b>x</b>	<b>x</b>		<b>x</b>	<b>x</b>	
<b>AB5019</b>	Sponge pads for tool AB5012, set 5 pcs								<b>o</b>
<b>VF2316</b>	TQC opacity scrub charts, black and white, 100 pieces	<b>o</b>			<b>o</b>	<b>o</b>			
<b>VF2316</b>	TQC opacity scrub charts black, 100 pieces.								

**Spare parts**

<b>AB5140</b>	Glass container 500ml								
<b>AB5142</b>	Glass container cap yellow								
<b>AB5141</b>	Glass container cap green								
<b>AB5113</b>	Pump tubes set of 4								
<b>AB5113</b>	Full tubing kit								